MUZAFFARPUR INSTITUTE OF TECHNOLOGY

Muzaffarpur

LAB MANUAL

Branch: Civil Engineering Year & Semester: 2nd Year/ 3rd Sem

FLUID MECHANICS LABORATORY

(**01 1X07**)



Prepared By:

Niraj Kumar Assistant Professor Department of Civil Engineering Atul Kumar Rahul Assistant Professor Department of Civil Engineering

FLUID MECHANICS LABORATORY

•

(01 1X07 P)

INDEX

Sl. No.	Name of the experiment	Date	Page no	Remarks
01	Stability of floating bodies			
02	Flow through notches			
03	To verify Bernoulli's theorem			
04	Impact of jet			
05	Coefficient of discharge, contraction & velocity of an Orifice			
06	Determination of Reynolds number			

Title: Stability of Floating Bodies

Aim:- To determine the Meta-centric height of a floating body.

Apparatus Used:- Take tank 2/3 full of water, floating vessel or pontoon fitted with a pointed pointer moving on a graduated scale, with weights adjusted on a horizontal beam.

Theory: -



Consider a floating body which is partially immersed in the liquid, when such a body is tilted, the center of buoyancy shifts from its original position 'B' to 'B' (The point of application of buoyanant force or upward force is known as center of G which may be below or above the center of buoyancy remain same and couple acts on the body. Due to this couple the body remains stable. At rest both the points G and B also Fb x Wc act through the same vertical line but in opposite direction. For small change (θ) B shifted to B.

The point of intersection M of original vertical line through B and G with the new vertical, line passing through 'B' is known as metacentre. The dis tance between G and M is known as metacentre height which is measure of static stability.

Date:

Formula Used:-

$$GM = \frac{W_m \cdot X_d}{(W_c + W_m) \tan \theta}$$

Where: -

•

Wm is unbalanced mass or weight.

Wc is weight of pontoon or anybody.

Xd is the distance from the center of pointer to striper or unbalanced weight. θ is angle of tilt or heel.

Procedure: -

- 1. Note down the dimensions of the collecting tank, mass density of water.
- 2. Note down the water level when pontoon is outside the tank.
- 3. Note down the water level when pontoon is inside the tank and their difference.
- 4. Fix the strips at equal distance from the center.
- 5. Put the weight on one of the hanger which gives the unbalanced mass.
- 6. Take the reading of the distance from center and angle made by pointer on arc.
- 7. The procedure can be repeated for other positioned and values of unbalanced mass.

Observation Table:-

Length of the tank = Width of the tank = Area of the tank = Initial level of the water without pontoon X1 = Final level of the water with pontoon (after adding unbalanced weight) X2 = Difference in height of water (X) = X2-X1=

Height of						
water in	Difference	Weight of	Tubalanaad		Analast	CM Mata acentria
tank with	in height	pontoon	Undalanced	Xd (m)	Angle of	GM=Metacentric
pontoon	X=X2-X1	Wc=XAp	masswm Kg		turn(Θ)	Height (m)
X2						

Precautions: -

•

- 1. The reading taking carefully without parallax error.
- 2. Put the weight on the hanger one by one.
- 3. Wait for pontoon to be stable before taking readings.
- 4. Strips should be placed at equal distance from the centre.

Result:- Meta centric height of the pontoon is measured with different positions and weights and value is.....

Viva Questions:-

- 1. Define Buoyancy?
- 2. Define Meta-centre?
- 3. Define Meta- centric height?
- 4. With respect to the position of metacentre, state the condition of equilibrium for a floating body?

•

Title: Flow through notches

Aim: - To determine the coefficient of discharge of V- Notch

Apparatus Used:- Arrangement for finding the coefficient of discharge inclusive of supply tank, collecting tank, pointer, scale & different type of notches

Theory:- Notches are overflow structure where length of crest along the flow of water is accurately shaped to calculate discharge.

Formula Used:-

For V notch the discharge:

Actual discharge Qa = $AY/t m^3/s$

Theoretical discharge $Q_t = 8/15\sqrt{2g}$. $\tan\theta/2 H^{5/2} m^3/s$

Co efficient of discharge $C_d = Q_{actual} / Q_{theoretical}$

Where:-

Q = DischargeH =Height above crest level $\theta = Angle of notch$ B = Width of notch

Procedure:-

`

- 1. The notch under test is positioned at the end of tank with vertical sharp edge on the upstream side.
- 2. Open the inlet valve and fill water until the crest of notch.
- 3. Note down the height of crest level by pointer gauge.
- 4. Change the inlet supply and note the height of this level in the tank.
- 5. Note the volume of water collected in collecting tank for a particular time and find out the discharge.
- 6. Height and discharge readings for different flow rate are noted.

Observations:-

Breath of tank = Length of tank = Height of water to crest level for V notch = Angle of V notch =

	Ι	Discharge Cal	lculation				Head		
Initial height of tank (Y ₁)	Final height of tank (Y ₂)	Difference in height (Y ₂ - Y ₁)	Volume (V)	Time (t)	Q = (V/t)	Height of water to crest level for V notch (H1)	Final height reading above crest (H ₂)	Head above crest level H= (H ₂₋ H ₁)	Cd

Precaution:-

•

- 1. Make the water level surface still, before takings the reading.
- 2. Reading noted should be free from parallax error.
- 3. The time of discharge is noted carefully.
- 4. Only the internal dimensions of collecting tank should be taken for consideration and calculations.

Result:

The value of C_d for V-Notch.....

Viva Questions:-

1. Differentiate between :-Uniform and non-

uniform flow Steady

and unsteady flow

- 2. Define notch?
- 3. What is coefficient of discharge?

Title: Verify the Bernoulli's Theorem

Aim:- To verify the Bernoulli's theorem.

Apparatus Used:- A supply tank of water, a tapered inclined pipe fitted with no. of piezometer tubes point, measuring tank, scale, stop watch.

Theory:- Bernoulli's theorem states that when there is a continues connection between the particle of flowing mass liquid, the total energy of any sector of flow will remain same provided there is no reduction or addition at any point.



Date:

Formula Used:-

`

$$H = \frac{P}{\rho g} + \frac{V^2}{2g} + Z$$

H = Total Head

Z= Potential head

$$P/\rho g = Pressure head$$

- $V^2/2g =$ Velocity head
- P = intensity of pressure

S = density of water

- g = acceleration due to gravity
- V = velocity of flowing water

$$H1 = H2 = H3 = H$$
 (Total Head)

Procedure:-

- 1. Open the inlet valve slowly and allow the water to flow from the supply tank.
- 2. Now adjust the flow to get a constant head in the supply tank to make flow in and out flow equal.
- 3. Under this condition the pressure head will become constant in the piezometer tubes.
- 4. Note down the quantity of water collected in the measuring tank for a given interval of time.
- 5. Compute the area of cross-section under the piezometer tube.
- 6. Compute the area of cross-section under the tube.
- 7. Change the inlet and outlet supply and note the reading.
- 8. Take at least three readings as described in the above steps.

Observation table:-

•

	1	2	3	4	5	6	7	8	9	10	11
Reading of piezometric tubes											
Area of cross section under the foot of each point											
Velocity of water under foot of each point											
$V^2/2g$											
p/ρ $V^2/2\alpha + p/\alpha$											
$v/2g + p/\rho$											

Sample calculation:

Dimensions of measuring tank

L = and B =

Area of measuring tank, (A) = $(L \times B)$ = Rise in water level in the measuring tank, (x)= Time taken for the water rise in the tank (t) = Actual Discharge (Q) = Area of cross section of the duct (a) = velocity (v) = Q/a velocity head = V²/2g (g = 9.81m/s)

Total Head (H) =
$$\frac{P}{\rho g} + \frac{V^2}{2g} + Z$$

Precautions:-

•

- 1. When fluid is flowing, there is a fluctuation in he height of piezometer tubes, note the mean position carefully.
- 2. Carefully keep some level of fluid in inlet and outlet supply tank.

Result:-

Viva Questions:-

- 1. Briefly explain the various terms involved in Bernoulli's equation?
- 2. Assumption made to get Bernoulli's equation from Euler's equation by made?

Date:

Title: Impact of jet

Aim: - To determine the co efficient of impact for vanes.

Apparatus Used:- Collecting tank, Transparent cylinder, Two nozzles of dia 10 mm & 12mm, Vane of different shape (flat, inclined or curved)

Theory:- Momentum equation is based on Newton's second law of motion which states that the algebraic sum of external forces applied to control volume of fluid in any direction is equal to the rate of change of momentum in that direction. The external forces include the component of the weight of the fluid & of the forces exerted externally upon the boundary surface of the control volume. If a vertical water jet moving with velocity is made to strike a target, which is free to move in the vertical direction then a force will be exerted on the target by the impact of jet, according to momentum equation this force (which is also equal to the force required to bring back the target in its original position) must be equal to the rate of change of momentum of the jet flow in that direction.



Figure: Impact of jet

Formula Used:-

`

 $F = \rho Q V(1 - \cos\beta)$ $F' = \rho Q 2 (1 - \cos\beta)$

 β =angle of difference vane

 ρ = density of water

Where F' =force (calculated)

V =velocity of jet angle

Q = discharge

A = area of nozzle ($\pi/4d2$)

For flat vane:

 $B = 90^{\circ}$

 $F = \rho Q^2 / a$

Procedure:-

- 1. Note down the relevant dimension or area of collecting tank, dia of nozzle, and density of water.
- 2. Install any type of vane i.e. flat, inclined or curved.
- 3. Install any size of nozzle i.e. 10mm or 12mm dia.
- 4. Note down the position of upper disk, when jet is not running.
- 5. Note down the reading of height of water in the collecting tank.
- 6. As the jet strike the vane, position of upper disk is changed, note the reading in the scale to which vane is raised.
- 7. Put the weight of various values one by one to bring the vane to its initial position.
- 8. At this position finds out the discharge also.
- 9. The procedure is repeated for each value of flowrate by reducing the water supply.

10. This procedure can be repeated for different type of vanes and nozzle.

Observation table:-

Dia of nozzle =

`

Mass density of water ρ =

Area of collecting tank =

Area of nozzle =

Horizontal flat vane:

When jet is not running, position of upper disk is at=

					Balancin				
Sr.	Discharg	e measure	ement		g		Theoretical	Error	in%
No.	Initial	Final	Time	Discharge	Mass	Force	force	=(F-F')	/F'
	(cm.)	(cm.)	(sec)	Q (cm ³ /sec)	(gm)	F(dyne)	$F' = \rho Q^2 / a$		
1.									
2.									
3.									

Result:

Precautions:-

- 1. Water flow should be steady and uniform.
- 2. The reading on the scale should be taken withoutany error.
- 3. The weight should be put slowly & one by one.
- 4. After changing the vane the flask should be closed tightly.

Viva Questions:-

- 1. Define the terms impact of jet and jet propulsion?
- 2. Find the expression for efficiency of a series of moving curved vane when a jet of water strikes the vanes at one of its tips?

Date:

Title: Coefficient of discharge, contraction & velocity of an Orifice

Aim:- To determine the coefficient of discharge, contraction & velocity of an Orifice.

Apparatus Used:- Supply tank with overflow arrangement, Orifice plate of different diameter, hook gauge, collecting tank, piezometric tube.

Theory:- A mouthpiece is a short length of pipe which is two or three times its diameter in length. If there pipe is filled externally to the orifices, the mouthpiece is called external cylindrical mouthpiece and discharge through orifice increase is a small opening of any cross-section on the side of bottom of the tank, through which the fluid is flowing orifice coefficient of velocity is defined as the ratio of two actual discharge to orifice ratio of the actual velocity of the jet at vena- contracta to the coefficient of theoretical velocity of the jet coefficient of contraction of defined as ratio of the actual velocity of jet at vena- contracta.

Vena- Contracta:- The fluid out is in form of jet goes on contracting form orifice up to dispute of about ¹/₂ the orifice dia. After the expend this least relation

Coefficient of velocity:- It is a ratio of actual velocity jet at vena-contracta to theoretical velocity.

Formula Used:-

$$C_d = \frac{Q_{actual}}{Q_{theoretical}}$$

 $Q_{\text{theoretical}} = \text{Theoretical velocity x Theoretical area} = \sqrt{2gh} \cdot a$

$$C_d = \frac{Q}{a \cdot \sqrt{2gh}}$$

 $C_v = rac{Actual \ velocity \ of \ jet \ at \ Vena \ contracta}{Theoretical \ Velocity}$

 $Coefficient of Contraction = \frac{Area of jet at Vena contracta}{Area of Orifice}$

$$C_c = \frac{a_c}{a}$$



Fig: Flow through Orifice



Procedure:-

•

- 1. Set the mouthpiece of orifice of which the Cc, Cu, Cd are to be determined.
- 2. Note the initial height of water in the steady flow tank and the height of datum from the bottom of orifice and mouthpiece. These remains constant for a particular mouthpiece or orifice.
- 3. By using the stop valve, set a particular flow in tank and tank height of water in tank.
- 4. Take the reading of discharge on this particular flow.
- 5. Using hook gauge, find the volume of Xo Y for mouthpiece.
- 6. Take three readings using hook gauge for one particular orifice.
- 7. Using the formula get value of Cd, Cu, and Cc for a particular orifice and mouthpiece.

Observation:-

Area of measuring tank =

Dia. Of orifice d =

Area of Orifice =

Discharge Calculation:

Sl. No.	Initial height h1	Final height h2	h=h2-h1	Time (t)	Discharge(Q=Ah/t)	

Co-efficient Calculation:

	Х	- axis			Y- ax	is				
Head			X=			V-V2	Qth =	Cd=	Cv=	Cc=
(H)	X1	X2	X2-	Y1	Y2	1-12 V1	$A\sqrt{2gh}$	Qact/Qth	$\sqrt{X^2/4}$ YH	Cd/Cv
			X1			-11				

Result:-

•

Precautions:-

- 1. Take the reading of discharge accurately.
- 2. Take value of h without any parallax error.
- 3. Set the orifice and mouthpiece.
- 4. Height of water in the steady flow.
- 5. Take reading from hook gauge carefully

Date:

Title: Determination of Reynolds number

Aim : To determine the Reynolds no with hele shaw apparatus

Apparatus Used: hele shaw Instrument, Measuring Cylinder, KMnO₄ Solution, Cylinder Container

Theory: This experiment employs the Hele shaw setup which produces a Flow pattern similar to that of potential flow. The flow is actually a highly viscous flow between two parallel plates with a very small gap between them. The flow through this apparatus is 2-D, low-speed. Although the flow is at low Reynolds number, this has a wide application in the Flow visualization apparatus as it produces the streamlines of potential flow.



Fig: Hele Shaw Apparatus

Formula Used:-

•

$$Q = \frac{AH}{t}$$

Q= Discharge

A= Cross-section Area of Measuring Cylinder

H= Difference in height

$$V = \frac{Q}{A}$$

V= Velocity

Q= Discharge

A= Area of the flow passage between the parallel plates from the tank (Area of slit)

Reynolds number (Re) =
$$\frac{\rho VL}{\mu}$$

 ρ = Mass density of fluid

V = Velocity

L = Characteristic Linear dimension

 μ = Dynamic viscosity

Procedure:

- 1. Mount the Hele-Shaw apparatus in a vertical position.
- 2. Place the model at the middle of the test section
- 3. Start the flow around the model
- 4. Time taken to fill the height of the measuring beaker shall be noted to calculate the volumetric flow rate
- 5. Photograph the flow

- 6. Change the angle of attack and visualize the change in flow features
- 7. Visualize the flow patterns for different models
- 8. Calculate Reynolds number

•

Observation table for Circular plate

Sl. No.	Height	Time(t)	Discharge (m ³ /s)	Area (m ²)	Velocity(m/s)	Chord Length (m)	Re
1.							
2.							
3.							

Observation table for Rectangular plate

Sl. No.	Height	Time(t)	Discharge (m ³ /s)	Area (m ²)	Velocity(m/s)	Chord Length	Re
						(m)	
1.							
2.							
3.							

Result: Reynolds number for rectangular and circular plate is------